


## User Guide to View Medical Images (USB)

The USB thumb drive issued by NUH DDI contains DICOM-quality images of radiology studies, as well as an integrated viewer software. Please follow these steps to view the study:

### Step 1: Insert the USB Thumb Drive

- Place the thumb drive into the computer USB port.
- Select decrypt\_win for Windows PC or decrypt\_mac for MacOS

Name	Date modified	Type	Size
decrypt_mac.app	20/07/2023 05:48 PM	File folder	
util	20/07/2023 05:48 PM	File folder	
VIEWER	20/07/2023 05:48 PM	File folder	
autorun.inf	20/07/2023 05:48 PM	Setup Information	1 KB
burn.enc	20/07/2023 05:47 PM	ENC File	77,722 KB
burn.lock	20/07/2023 05:47 PM	LOCK File	0 KB
 decrypt_win.exe	07/05/2023 02:46 PM	Application	11,839 KB
run.bat	31/07/2022 10:00 PM	Windows Batch File	1 KB

### Step 2 : Key in password when prompted.

- The password format will be:
  - 1)"DDI"
  - 2)"YY"**(Year of scan)**
  - 3)"Initials"(Follow the printed label, Jeremy Tan Wei Jie = JTWJ)
  - 4)"!!"

Example : DDI26JTWJ!!

